

## SPECIFICATION SHEET


### PRODUCT NAME:

ETFE Ferrule (natural) with SS lock ring for 1/16" (1.6 mm) OD Tubing



\*Image for reference only. Color of ferrule will vary as indicated in the product name.

### DIMENSIONS & TOLERANCE:



Tolerance is  $\pm 0.1$ mm unless otherwise specified. SS Ring dimensions not available.

### MAX OPERATING CONDITIONS:

Temperature: 70 °C

Note: Max temperature rating may vary depending on liquid.

Pressure Rating: 9 MPa (Liquid)\*; 900 kPa (Air)

\*Tested with water at 25 °C

### PACKAGING:

Distributed in bags with appropriate identifying label. Custom packaging and labeling are available upon request.

### COMPLIANCE INFORMATION:

Compliant with RoHS and REACH regulation.

### REVISION HISTORY

REV	Release Date	Description
A	11/05/25	Initial Release.

END OF DOCUMENT