

SPECIFICATION SHEET


PRODUCT NAME:

ETFE Ferrule (natural) with SS lock ring for 1/8" (3.2 mm) OD Tubing



*Image for reference only. Color of ferrule will vary as indicated in the product name.

DIMENSIONS & TOLERANCE:



Tolerance is $\pm 0.1\text{mm}$ unless otherwise specified. SS Ring dimensions not available.

MAX OPERATING CONDITIONS:

Temperature: 70 °C

Note: Max temperature rating may vary depending on liquid.

Pressure Rating: 6 MPa (Liquid)*; 700 kPa (Air)

*Tested with water at 25 °C

PACKAGING:

Distributed in bags with appropriate identifying label. Custom packaging and labeling are available upon request.

COMPLIANCE INFORMATION:

Compliant with RoHS and REACH regulation.

REVISION HISTORY

REV	Release Date	Description
A	11/05/25	Initial Release.

END OF DOCUMENT